

Rapid Barcoding Kit 96 (SQK-RBK110.96)

Version: RBK_9126_v110_revB_24Mar2021
Last update: 21/04/2021



Flow Cell Number:

DNA Samples:

Before start checklist

Materials

- ☐ 50 ng high molecular weight genomic DNA per sample
- ☐ For R10.3 flow cells, 50-200 ng high molecular weight genomic DNA per sample
- ☐ Rapid Barcoding Kit 96 (SQK-RBK110.96)

Consumables

- ☐ 1.5 ml Eppendorf DNA LoBind tubes
- ☐ 2 ml Eppendorf DNA LoBind tubes
- ☐ 0.2 ml thin-walled PCR tubes
- ☐ Eppendorf twin.tec® PCR plate 96 LoBind, semi-skirted (Cat # 0030129504) with heat seals
- ☐ Nuclease-free water (e.g. ThermoFisher, cat # AM9937)
- ☐ Freshly-prepared 80% ethanol in nuclease-free water

Equipment

- ☐ Ice bucket with ice
- ☐ Microplate centrifuge, e.g. Fisherbrand™ Mini Plate Spinner Centrifuge (Fisher Scientific, #11766427)
- ☐ Timer
- ☐ Thermal cycler or heat block at 30°C and 80°C
- ☐ Magnetic rack
- ☐ Hula mixer (gentle rotator mixer)
- ☐ Pipettes and pipette tips P2, P20, P100, P200, P1000

INSTRUCTIONS

NOTES/OBSERVATIONS

Library preparation

- ☐ Program the thermal cycler: 30°C for 2 minutes, then 80°C for 2 minutes.
- Thaw kit components at RT, spin down briefly using a microfuge and mix by pipetting as indicated by the table below:
- ☐ Rapid Barcode plate (RB96): not frozen, briefly spin down, mix well by pipetting
 - ☐ Rapid Adapter F (RAP-F): not frozen, briefly spin down, mix well by pipetting
 - ☐ SPRI beads (SPRI): thaw at RT, briefly spin down, mix by pipetting or vortexing immediately before use
 - ☐ Sequencing Buffer II (SBII): thaw at RT, briefly spin down, mix well by pipetting*
 - ☐ Loading Beads II (LBII): thaw at RT, briefly spin down, mix by pipetting or vortexing immediately before use
 - ☐ Elution Buffer (EB): thaw at RT, briefly spin down, mix well by pipetting
 - ☐ Flush Buffer (FB): thaw at RT, briefly spin down, mix by vortexing
 - ☐ Flush Tether (FLT): thaw at RT, briefly spin down, mix well by pipetting
- Prepare the DNA in Nuclease-free water.
- ☐ Transfer 50 ng genomic DNA per sample (if using R9.4.1 flow cells) or 50-200 ng (if using R10.3 flow cells) into a 1.5 ml Eppendorf DNA LoBind tube
 - ☐ Adjust the volume to 9 µl with Nuclease-free water
 - ☐ Mix by pipetting up and down
 - ☐ Spin down briefly in a microfuge

Rapid Barcoding Kit 96 (SQK-RBK110.96)

Version: RBK_9126_v110_revB_24Mar2021
Last update: 21/04/2021



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>In 0.2 ml thin-walled PCR tubes or an Eppendorf twin.tec® PCR plate 96 LoBind, mix the following. The Rapid Barcodes can be transferred using a multichannel pipette:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 9 µl 50 ng template DNA <input type="checkbox"/> 1 µl Rapid Barcodes (RB01-96, one for each sample) <p><input type="checkbox"/> Mix well by pipetting. Seal the plate and spin down in a centrifuge.</p> <p><input type="checkbox"/> Incubate the tubes or plate at 30°C for 2 minutes and then at 80°C for 2 minutes. Briefly put the tubes or plate on ice to cool.</p> <p><input type="checkbox"/> Pool all barcoded samples in your desired ratio, noting the total volume.</p> <p><input type="checkbox"/> Resuspend the SPRI beads by vortexing.</p> <p><input type="checkbox"/> To the entire pooled barcoded sample from Step 7, add an equal volume of resuspended SPRI beads and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare at least 3 ml of fresh 80% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 1.5 ml of freshly-prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Briefly spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for 30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 15 µl Elution Buffer (EB). Incubate for 10 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p>Remove and retain 15 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Dispose of the pelleted beads <p><input type="checkbox"/> Quantify 1 µl of eluted sample using a Qubit fluorometer.</p> <p><input type="checkbox"/> Transfer 11 µl of the sample into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Add 1 µl of Rapid Adapter F (RAP-F) to 11 µl of barcoded DNA.</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the reaction for 5 minutes at RT.</p>	
<p>The prepared library is used for loading into the MinION flow cell. Store the library on ice until ready to load.</p>	

Rapid Barcoding Kit 96 (SQK-RBK110.96)

Version: RBK_9126_v110_revB_24Mar2021
Last update: 21/04/2021



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Priming and loading the SpotON Flow Cell</p> <p>Using the Loading Solution</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer II (SBII), Loading Beads II (LBII) or Loading Solution (LS, if using), Flush Tether (FLT) and Flush Buffer (FB) at RT. <input type="checkbox"/> Mix the Sequencing Buffer II (SBII), Flush Buffer (FB) and Flush Tether (FLT) tubes by vortexing. Spin down the SBII and FLT at RT. <input type="checkbox"/> Open the MinION Mk1B lid and slide the flow cell under the clip. <input type="checkbox"/> Slide the priming port cover clockwise to open the priming port. <p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. <p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles (a few µl):</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <ul style="list-style-type: none"> <input type="checkbox"/> To prepare the flow cell priming mix, add 30 µl of thawed and mixed Flush Tether (FLT) to 1.17 ml of thawed and mixed Flush Buffer (FB), and mix by vortexing at RT. <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes. During this time, prepare the library for loading by following the steps below. <input type="checkbox"/> Thoroughly mix the contents of the Loading Beads II (LBII) by pipetting. <p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The Loading Beads II (LBII) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use. <p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl Sequencing Buffer II (SBII) <input type="checkbox"/> 25.5 µl Loading Beads II (LBII) mixed immediately before use, or Loading Solution (LS), if using <input type="checkbox"/> 12 µl DNA library <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. 	

Rapid Barcoding Kit 96 (SQK-RBK110.96)

Version: RBK_9126_v110_revB_24Mar2021
Last update: 21/04/2021



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION Mk1B lid.	
Ending the experiment	
<input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Wash Kit instructions and store the washed flow cell at 2-8°C, OR	
<input type="checkbox"/> Follow the returns procedure by washing out the flow cell ready to send back to Oxford Nanopore.	
IMPORTANT	
<input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.	