

# Ultra-Long DNA Sequencing Kit (SQK-ULK001)

Version: ULK\_9124\_v110\_revA\_24Mar2021  
 Last update: 24/03/2021



Flow Cell Number: .....

DNA Samples: .....

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 750 µl of extracted UHMW gDNA in Circulomics EB+	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> Ultra-Long DNA Sequencing Kit (SQK-ULK001)		<input type="checkbox"/> Thermal cycler or heat block at 75°C
<input type="checkbox"/> Circulomics Nanobind Big DNA Kit (e.g. CBB, Plant or Tissue)		<input type="checkbox"/> Vortex mixer
<input type="checkbox"/> Circulomics Nanobind UL Library Prep Kit (NB-900-601-01)		<input type="checkbox"/> Microfuge
		<input type="checkbox"/> Ice bucket with ice
		<input type="checkbox"/> Timer
		<input type="checkbox"/> Pipettes and pipette tips P200

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Tagmentation</b></p> <p>Before starting the library preparation, extract 750 µl of UHMW gDNA in Circulomics EB+.</p> <p><input type="checkbox"/> Thaw the Fragmentation Mix (FRA), FRA Dilution Buffer (FDB) and Rapid Adapter-F (RAP-F). Spin down briefly using a microfuge and keep on ice.</p> <p>In a 1.5 ml Eppendorf DNA LoBind tube, make up the diluted FRA:</p> <p><input type="checkbox"/> 6 µl Fragmentation Mix (FRA)</p> <p><input type="checkbox"/> 244 µl FRA dilution buffer (FDB)</p> <p><input type="checkbox"/> Mix the diluted FRA by vortexing.</p> <p><input type="checkbox"/> Add 250 µl of diluted FRA to the extracted DNA. Stir the reaction with the pipette tip whilst expelling the diluted FRA to ensure an even distribution.</p> <p><input type="checkbox"/> Immediately vortex the reaction for 5 seconds at the lowest setting able to generate a gentle vortex in the fluid.</p> <p>Incubate the reaction as follows:</p> <p><input type="checkbox"/> 5 minutes Room temperature</p> <p><input type="checkbox"/> 5 minutes 75°C</p> <p><input type="checkbox"/> Cool to RT for a minimum of 10 minutes Room temperature</p> <p><input type="checkbox"/> Add 5 µl of RAP-F with a regular pipette tip. Use a P1000 wide-bore tip to pipette mix. Visually check to ensure the reaction is thoroughly mixed. Inversion can be used to aid mixing.</p> <p><input type="checkbox"/> Incubate for 30 minutes at RT.</p>	

# Ultra-Long DNA Sequencing Kit (SQK-ULK001)

Version: ULK\_9124\_v110\_revA\_24Mar2021  
 Last update: 24/03/2021



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Clean-up</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Thaw one tube of Long Fragment Buffer (LFB) and Elution Buffer (EB) at RT, mix by vortexing, spin down and store at RT.</li> <li><input type="checkbox"/> Add a single Circulomics Nanobind disk to the DNA library and invert to ensure both sides of the disk are wetted.</li> <li><input type="checkbox"/> Add 500 µl Circulomics NAF Binding Buffer.</li> <li><input type="checkbox"/> Gently invert the tube 20-30 times to mix until all reagents are thoroughly mixed.</li> <li><input type="checkbox"/> Wait 5-10 minutes. The DNA should start to become visible as a precipitate forming around the disk.</li> <li><input type="checkbox"/> Invert the tube again 5-10 times to ensure all DNA has precipitated and is tightly bound to the Nanobind Disk.</li> <li><input type="checkbox"/> Place the tube on a magnetic rack to aspirate and discard the supernatant with a P200 pipette. Take care not to disturb the DNA precipitated onto the Nanobind disk.</li> <li><input type="checkbox"/> Wash the Nanobind disk by gently adding 1000 µl of Long Fragment Buffer (LFB) and invert. Keep the tube on the magnetic rack and incubate for 5 minutes at RT. Aspirate and discard the supernatant, taking care not to disturb the DNA precipitate.</li> <li><input type="checkbox"/> Repeat the previous step and remove any supernatant in the lid of the Eppendorf tube.</li> <li><input type="checkbox"/> Add 225 µl of Elution Buffer (EB) to elute the DNA.</li> <li><input type="checkbox"/> Incubate the library overnight at RT, for a minimum of 12 hours. Gently aspirate and dispense the eluate over the Nanobind disk at regular intervals with a wide-bore pipette tip to aid elution.</li> <li><input type="checkbox"/> Transfer 225 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube after the overnight incubation with a wide-bore tip. The DNA should be viscous.</li> <li><input type="checkbox"/> Spin the tube containing the Nanobind disk on a centrifuge at 10,000 x g for 10-15 seconds and transfer any additional liquid that comes off the disk to the previous eluate. Repeat as necessary until all DNA is removed using a standard P200 pipette.</li> <li><input type="checkbox"/> Gently mix the eluate containing the DNA library five times with a wide-bore pipette tip.</li> <li><input type="checkbox"/> Incubate for 2 hours at RT. The library should be very viscous.</li> </ul>	
<p>The prepared library is used for loading into the flow cell. Store at RT to allow DNA to homogenise before loading.</p>	
<p><b>Priming and loading the SpotON flow cell for GridION</b></p>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Priming Tether (PT) and one tube of Flush Buffer (FB) at RT. Mix by vortexing and spin down.</li> </ul> <p>In a new tube, prepare the DNA library for loading as follows and gently mix with a wide-bore tip:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 37.5 µl Sequencing Buffer (SQB)</li> <li><input type="checkbox"/> 37.5 µl DNA library</li> </ul>	

# Ultra-Long DNA Sequencing Kit (SQK-ULK001)

Version: ULK\_9124\_v110\_revA\_24Mar2021  
 Last update: 24/03/2021



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <li><input type="checkbox"/> Wait 30 minutes and gently mix with a wide-bore tip. Visually inspect to ensure the sample is homogenous.</li> <li><input type="checkbox"/> Slide open the GridION lid and insert the flow cell.</li> <li><input type="checkbox"/> Slide the priming port cover clockwise to open the priming port.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 <math>\mu</math>l, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</li> </ul>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles (a few <math>\mu</math>l):</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 <math>\mu</math>l</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 <math>\mu</math>l, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Prepare the flow cell priming mix in a 1.5 ml Eppendorf tube and mix by vortexing at RT.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 30 <math>\mu</math>l Priming Tether (PT)</li> <li><input type="checkbox"/> 1170 <math>\mu</math>l Flush Buffer (FB)</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Load 800 <math>\mu</math>l of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes.</li> </ul> <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li> <li><input type="checkbox"/> Load 200 <math>\mu</math>l of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> <li><input type="checkbox"/> Thoroughly mix the DNA library and ensure the SpotON port and Priming port covers of the flow cell are open in preparation for loading.</li> <li><input type="checkbox"/> Load the DNA library dropwise onto the SpotON port until 75 <math>\mu</math>l has been loaded. Ensure each drop flows into the port before adding the next.</li> <li><input type="checkbox"/> Cover Waste port 2 and the Priming port with clean, gloved fingers.</li> <li><input type="checkbox"/> Using a fully depressed P200 pipette, insert the tip in Waste port 1 whilst Waste port 2 and the Priming port are covered.</li> <li><input type="checkbox"/> Very slowly aspirate to pull the DNA library into the SpotON sample port. Closely watch the DNA library on the SpotON port and completely remove the pipette as soon as the library starts to be pulled into the port.</li> <li><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the GridION lid.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> When setting up your run parameters in MinKNOW, under Run Options, set the time between mux scans to 6 hours from the default 1.5 hours before starting the sequencing run.</li> </ul>	
<p>We recommend loading an ultra-long DNA library three times per flow cell to increase output.</p>	

# Ultra-Long DNA Sequencing Kit (SQK-ULK001)

Version: ULK\_9124\_v110\_revA\_24Mar2021  
 Last update: 24/03/2021



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Reloading ultra-long DNA library on a MinION/GridION flow cell</b></p>	
<p>To run a second library of ultra-long DNA straight after flushing a flow cell, we recommend removing all fluid from the waste channel after each priming step, as outlined below.</p> <p><input type="checkbox"/> Mix the Sequencing Buffer (SQB), Priming Tether (PT) and Flush Buffer (FB) tubes by vortexing and spin down at RT.</p> <p>In a new tube, prepare the DNA library for loading as follows and gently mix with a wide-bore tip:</p> <p><input type="checkbox"/> 37.5 µl Sequencing Buffer (SQB)</p> <p><input type="checkbox"/> 37.5 µl DNA library</p> <p><input type="checkbox"/> Wait 30 minutes and gently mix with a wide-bore tip. Visually inspect to ensure the sample is homogenous.</p> <p><input type="checkbox"/> Slide the priming port cover of the flow cell clockwise to open the priming port.</p>	
<p><b>IMPORTANT</b></p> <p><input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.</p>	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles (a few µl):</p> <p><input type="checkbox"/> Set a P1000 pipette to 200 µl</p> <p><input type="checkbox"/> Insert the tip into the priming port</p> <p><input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip</p> <p>Prepare the flow cell priming mix in a 1.5 ml Eppendorf tube and mix by vortexing at RT.</p> <p><input type="checkbox"/> 30 µl Priming Tether (PT)</p> <p><input type="checkbox"/> 1170 µl Flush Buffer (FB)</p> <p><input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes.</p> <p><input type="checkbox"/> Close the priming port cover and ensure the SpotON sample port cover is closed.</p> <p><input type="checkbox"/> Using a P1000, remove all fluid from the waste channel through Waste port 1. As both the priming port and SpotON sample port are closed, no fluid should leave the sensor array area.</p> <p><input type="checkbox"/> Slide open the priming port and load 200 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles.</p> <p><input type="checkbox"/> Close the priming port and use a P1000 to remove all fluid from the waste channel through Waste Port 1.</p> <p><input type="checkbox"/> Thoroughly mix the DNA library and ensure the SpotON port and Priming port covers of the flow cell are open in preparation for loading.</p> <p><input type="checkbox"/> Load the DNA library dropwise onto the SpotON port until 75 µl has been loaded. Ensure each drop flows into the SpotON port before adding the next.</p> <p><input type="checkbox"/> Cover Waste port 2 and the Priming port with clean, gloved fingers.</p>	

# Ultra-Long DNA Sequencing Kit (SQK-ULK001)

Version: ULK\_9124\_v110\_revA\_24Mar2021  
Last update: 24/03/2021



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"><li><input type="checkbox"/> Using a fully depressed P200 pipette, insert the tip in Waste port 1 whilst Waste port 2 and the Priming port are covered.</li><li><input type="checkbox"/> Very slowly aspirate to pull the DNA library into the SpotON sample port. Closely watch the DNA library on the SportON port and completely remove the pipette as soon as the library starts to be pulled into the port.</li><li><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION Mk1B lid.</li><li><input type="checkbox"/> Once the flow cell is reloaded, resume the sequencing run on MinKNOW and trigger a mux scan.</li></ul>	
<b>Ending the experiment</b>	
<ul style="list-style-type: none"><li><input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Wash Kit instructions and store the washed flow cell at 2-8°C, OR</li><li><input type="checkbox"/> Follow the returns procedure by washing out the flow cell ready to send back to Oxford Nanopore.</li></ul>	
<b>IMPORTANT</b> <ul style="list-style-type: none"><li><input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.</li></ul>	